


# MIA – Image Formats

# Why DICOM?


2



# DICOM & PACS


3

- DICOM: Digital Imaging and COmmunications in Medicine
- PACS: Picture Archiving and Communication System



# DICOM & PACS

4



# Origins

5



NEMA, Suite 1752  
1300 North 17<sup>th</sup> Street  
Rosslyn, VA 22209  
Ph: (703) 841-3285  
<http://dicom.nema.org>

- ACR:
  - ▣ American College of Radiology
- NEMA:
  - ▣ National Electrical Manufacturers Association
- AAPM:
  - ▣ American Association of Physicists in Medicine
- RSNA:
  - ▣ Radiological Society of North America

# DICOM

6

- Proposed in 1983.
- Supports many imaging modalities

## CONTENTS

NOTICE AND DISCLAIMER.....	2
CONTENTS .....	3
FOREWORD.....	4
INTRODUCTION.....	5
HISTORY .....	5
THE DICOM STANDARD.....	5
CURRENT DIRECTION .....	6
RETIREMENT .....	6
1 Scope and field of application.....	7
2 Normative references .....	7
3 Definitions.....	8
4 Symbols and abbreviations .....	9
5 Goals of the DICOM standard .....	9
6 Overview of the Content of the DICOM Standard .....	11
6.1 DOCUMENT STRUCTURE.....	11
6.2 PS 3.2: CONFORMANCE .....	11
6.3 PS 3.3: INFORMATION OBJECT DEFINITIONS .....	14
6.4 PS 3.4: SERVICE CLASS SPECIFICATIONS.....	15
6.5 PS 3.5: DATA STRUCTURE AND SEMANTICS .....	15
6.6 PS 3.6: DATA DICTIONARY .....	16
6.7 PS 3.7: MESSAGE EXCHANGE.....	16
6.8 PS 3.8: NETWORK COMMUNICATION SUPPORT FOR MESSAGE EXCHANGE.....	16
6.9 PS 3.9: RETIRED (FORMERLY POINT-TO-POINT COMMUNICATION SUPPORT FOR MESSAGE EXCHANGE) .....	17
6.10 PS 3.10 MEDIA STORAGE AND FILE FORMAT .....	17
6.11 PS 3.11: MEDIA STORAGE APPLICATION PROFILES.....	18
6.12 PS 3.12: STORAGE FUNCTIONS AND MEDIA FORMATS FOR DATA INTERCHANGE.....	19
6.13 PS 3.13: RETIRED (FORMERLY PRINT MANAGEMENT POINT-TO-POINT COMMUNICATION SUPPORT).....	19
6.14 PS 3.14: GRAYSCALE STANDARD DISPLAY FUNCTION.....	20
6.15 PS 3.15: SECURITY AND SYSTEM MANAGEMENT PROFILES .....	20
6.16 PS 3.16: CONTENT MAPPING RESOURCE .....	20
6.17 PS 3.17: EXPLANATORY INFORMATION.....	20
6.18 PS 3.18: WEB ACCESS TO DICOM PERSISTENT OBJECTS (WADO).....	20
6.19 PS3.19: APPLICATION HOSTING .....	21
6.20 PS 3.20: TRANSFORMATION OF DICOM TO AND FROM HL7 STANDARDS .....	22



<http://medical.nema.org/standard.html>

# DICOM

7

- Information about
  - ▣ Patient: name, years, ...
  - ▣ Image acquisition: who, where, ...
  - ▣ Image information: resolution, spacing, ...
  - ▣ Image or Images
  - ▣ ...
  
- But also, the DICOM establishes a protocol related with the communication about the different DICOM applications (acquisition machines, PACS, reading stations, etc.)


## □ Hierarchical structure:





- The DICOM is structured in Data Elements, each of them numbered (Tag) and with a unique function:

Value Representation	Description
AE	Application Entity
AS	Age String
AT	Attribute Tag
CS	Code String
DA	Date
DS	Decimal String
DT	Date/Time
FL	Floating Point Single (4 bytes)
FD	Floating Point Double (8 bytes)
IS	Integer String
LO	Long String
LT	Long Text
OB	Other Byte
OF	Other Float
OW	Other Word
PN	Person Name
SH	Short String
SL	Signed Long
SQ	Sequence of Items
SS	Signed Short
ST	Short Text
TM	Time
UI	Unique Identifier
UL	Unsigned Long



# dicominfo('file.dcm')

10

**Filename:** [1x55char]  
**FileModDate:** '27-maig-200823: 25: 47'  
**FileSize:** 17042728  
**Format:** 'DICOM'  
**FormatVersion:** 3  
**Width:** 2560  
**Height:** 3328  
**BitDepth:** 12  
**ColorType:** 'grayscale'  
**FileMetaInformationGroupLength:** 212  
**FileMetaInformationVersion:** [2x1uint8]  
**MediaStorageSOPClassUID:** '1.2.840.10008.5.1.4.1.1.1.2'  
**MediaStorageSOPInstanceUID:** [1x55char]  
**TransferSyntaxUID:** '1.2.840.10008.1.2'  
**ImplementationClassUID:** '1.3.6.1.4.1.5962.99.2'  
**ImplementationVersionName:** 'PIXELMEDJAVA001'  
**SourceApplicationEntityTitle:** 'RXMAMOPACS'  
**SpecificCharacterSet:** 'ISO\_IR100'  
**ImageType:** 'DERIVED\PRIMARY\RIGHT'  
**SOPClassUID:** '1.2.840.10008.5.1.4.1.1.1.2'  
**SOPInstanceUID:** [1x55char]  
**StudyDate:** '20070803'  
**SeriesDate:** '20070803'  
**AcquisitionDate:** '20070803'  
**ContentDate:** '20070803'  
**StudyTime:** '120122.375000'  
**SeriesTime:** '120334.390000'  
**AcquisitionTime:** '120310.812000'

**ContentTime:** '120330.281000'  
**AccessionNumber:** '660471'  
**Modality:** 'MG'  
**PresentationIntentType:** 'FOR PRESENTATION'  
**Manufacturer:** 'SIEMENS'  
**InstitutionName:** 'HOSPITAL JOSEP TRUETA'  
**InstitutionAddress:** [1x27char]  
**ReferringPhysicianName:** [1x1struct]  
**StationName:** 'YBFR015811'  
**SeriesDescription:** 'Mamobilateral, Diagnose'  
**InstitutionalDepartmentName:** 'Department'  
**PerformingPhysicianName:** [1x2struct]  
**OperatorName:** [1x1struct]  
**ManufacturerModelName:** 'Mammomat Novation DR'  
**ReferencedPerformedProcedureStepSequence:** [1x1struct]  
**DerivationDescription:** [1x42char]  
**SourceImageSequence:** [1x1struct]  
**AnatomicRegionSequence:** [1x1struct]  
**PatientName:** [1x1struct]  
**PatientID:** '660471'  
**PatientBirthDate:** '19560506'  
**PatientSex:** 'F'  
**PatientAge:** '051Y'  
**ContrastBolusAgent:** ''  
**BodyPartExamined:** 'BREAST'  
**KVP:** 25  
**DeviceSerialNumber:** '1036'  
**SoftwareVersion:** [1x47char]

# dicominfo('file.dcm')... cont

11

## **DistanceSourceToDetector: 650**

DistanceSourceToPatient: 636

TableAngle: 0

TableType: 'NONE'

FieldOfViewShape: 'RECTANGLE'

FieldOfViewDimensions: [2x1double]

## **ExposureTime: 410**

XrayTubeCurrent: 147

Exposure: 60

ExposureInuAs: 60099

RectificationType: 'CONST POTENTIAL'

ImagerPixelSpacing: [2x1double]

Grid: 'FOCUSEDPARALLEL'

FocalSpot: 0.3000

## **AnodeTargetMaterial: 'TUNGSTEN'**

BodyPartThickness: 27

CompressionForce: 112

## **DateOfLastCalibration: '18590101'**

TimeOfLastCalibration: '000000.000000'

AcquisitionDeviceProcessingDescription: ''

AcquisitionDeviceProcessingCode: ''

TomoLayerHeight: 0

## **PositionerType: 'MAMMOGRAPHIC'**

PositionerPrimaryAngle: 0

PositionerSecondaryAngle: 0

DetectorPrimaryAngle: 0

DetectorSecondaryAngle: 0

ShutterShape: 'RECTANGULAR'

ShutterLeftVerticalEdge: 0

ShutterRightVerticalEdge: 2560

ShutterUpperHorizontalEdge: 0

ShutterLowerHorizontalEdge: 3328

ViewPosition: 'CC'

Sensitivity: 0

DetectorConditionsNominalFlag: 'YES'

DetectorTemperature: 30

## **DetectorType: 'DIRECT'**

DetectorConfiguration: 'AREA'

DetectorDescription: ''

DetectorMode: 'Normalmode'

DetectorID: 'MP1775'

DateOfLastDetectorCalibration: '20070724'

TimeOfLastDetectorCalibration: '082552.265000'

## **ExposuresOnDetectorSinceLastCalibration: 23702**

ExposuresOnDetectorSinceManufactured: 57839

DetectorTimeSinceLastExposure: 920.1090

DetectorActiveTime: 410

DetectorActivationOffsetFromExposure: 0

DetectorBinning: [2x1double]

DetectorElementPhysicalSize: [2x1double]

DetectorElementSpacing: [2x1double]

DetectorActiveShape: 'RECTANGLE'

DetectorActiveDimensions: [2x1double]

DetectorActiveOrigin: [2x1double]

FieldOfViewOrigin: [2x1double]

FieldOfViewRotation: 0

FieldOfViewHorizontalFlip: 'NO'

## **FilterMaterial: 'RHODIUM'**

# dicominfo('file.dcm')... cont



12

FilterThicknessMinimum: 0.0500  
FilterThicknessMaximum: 0.0500  
**ExposureControlMode: 'AUTOMATIC'**  
ExposureControlModeDescription: [1x55char]  
ExposureStatus: 'NORMAL'  
ExposureTimeInuS: 410000  
XrayTubeCurrentInuA: 1.4658e+005  
**StudyInstanceUID: [1x55char]**  
**SeriesInstanceUID: [1x55char]**  
StudyID: '1'  
SeriesNumber: 2  
AcquisitionNumber: 1  
InstanceNumber: 1  
PatientOrientation: 'P\L'  
**ImageLaterality: 'R'**  
ImageComments: ''  
SamplesPerPixel: 1  
PhotometricInterpretation: 'MONOCHROME2'  
Rows: 3328  
Columns: 2560  
BitsAllocated: 16  
BitsStored: 12  
HighBit: 11  
PixelRepresentation: 0  
QualityControlImage: 'NO'  
BurnedInAnnotation: 'NO'  
**PixelIntensityRelationship: 'LOG'**  
PixelIntensityRelationshipSign: 1

**WindowCenter: [2x1double]**  
**WindowWidth: [2x1double]**  
RescaleIntercept: 0  
RescaleSlope: 1  
RescaleType: 'US'  
ImplantPresent: 'NO'  
PartialView: 'NO'  
WindowCenterWidthExplanation: 'WINDOW1\WINDOW2'  
LossyImageCompression: '00'  
Private\_0029\_10xx\_Creator: 'SIEMENS MEDCOM HEADER'  
Private\_0029\_11xx\_Creator: 'SIEMENS MEDCOM HEADER2'  
**Private\_0029\_1031:** [12x1uint8]  
**Private\_0029\_1032:** [4x1uint8]  
**Private\_0029\_1033:** [4x1uint8]  
**Private\_0029\_1034:** [12x1uint8]  
Private\_0029\_1160: [4x1uint8]  
StudyPriorityID: 'MED'  
PerformedProcedureStepStartDate: '20070803'  
PerformedProcedureStepStartTime: '120122.375000'  
PerformedProcedureStepID: 'MG20070803120122'  
EntranceDose: 0  
DistanceSourceToEntrance: 609  
OrganDose: 0.0070  
**OrganExposed: 'BREAST'**  
AcquisitionContextSequence: [1x1struct]  
EntranceDoseInmGy: 2  
Private\_0041\_10xx\_Creator: 'SIEMENSMEDSPDXMGWHAWS1'  
Private\_0041\_1002: [10x1uint8]  
CalibrationImage: 'NO'  
ViewCodeSequence: [1x1struct]  
HistogramSequence: [1x1struct]  
Private\_0095\_10xx\_Creator: 'SIENET'  
Private\_0095\_100c: [4x1uint8]  
PresentationLUTShape: 'IDENTITY'



# dicomread('file.dcm')

13



# dicomread('file.dcm')

14



# ISSUES


15

- Ethics. Patient related data. Anonymize the data!
- Formats. 2D, 3D, 4D, vendor specific.
- Visualization, interpretation of the DICOM header
- Image bit depth (8, 12, 16)
- Image dimensionality
- Other formats [http://www.itk.org/Wiki/ITK/File\\_Formats](http://www.itk.org/Wiki/ITK/File_Formats)
  - ▣ nrrd (Nearly Raw raster data <http://teem.sourceforge.net/nrrd/>)
  - ▣ mhd (meta image)
  - ▣ Analyze, etc...

# Images & Pixels (from itk slides)

16

- What is an image?
  - ▣ Sampling of a continuous field using a discrete grid.






# Images & Pixels (from itk slides)

17


## □ Sampling grid



# Images & Pixels (from itk slides)

18

## □ Image indices



# Images & Pixels (from itk slides)

19

- Use pixel or physical coordinates?
- Convert from physical to pixel coordinates

$$P[0] = \text{Index}[0] \times \text{Spacing}[0] + \text{Origin}[0]$$

$$P[1] = \text{Index}[1] \times \text{Spacing}[1] + \text{Origin}[1]$$

$$\text{Index}[0] = \text{floor}( ( P[0] - \text{Origin}[0] ) / \text{Spacing}[0] + 0.5 )$$


$$\text{Index}[1] = \text{floor}( ( P[1] - \text{Origin}[1] ) / \text{Spacing}[1] + 0.5 )$$

- Activity: Physical coordinates example.
  - ▣ 3D Ultrasound. What is the size of the lesion?

# Images & Pixels (from itk slides)

20


## □ Image region



# Images & Pixels (from itk slides)

21


## □ Resampling



# Images & Pixels (from itk slides)

22


## □ Super Sampling



# Images & Pixels (from itk slides)

23

## □ Comparing images



# Things I will not do...

24

*I will not register images in pixel space*  
*I will not register images in pixel space*  
*I will not register images in pixel space*  
*I will not register images in pixel space*  
*I will not register images in pixel space*  
*I will not register images in pix*





# To know more...

25

- ❑ [medical.nema.org](http://medical.nema.org) – Official Site
- ❑ [rsna.org/Technology/DICOM](http://rsna.org/Technology/DICOM)
- ❑ [Wikipedia.org/DICOM](http://Wikipedia.org/DICOM)
- ❑ [Itk.org](http://Itk.org)
- ❑ Insight into images